

New instruments in Implantology:

THE PPBUSER SCRAPERS

PROF. DR. DANIEL BUSER

PPBUSER.....

These scrapers permit the use of autologous bones. The scrapers have rounded corners and are available in the following widths:

PPBUSE 1/26 3 + 4 mm

PPBUSE 1/36 3 + 5 mm

PPBUSE 2/36 4 + 5 mm




PPBUSE 1/26

With this scraper it is possible to remove excellent transplantation material in suitable bone quantities.

ADVANTAGES:

- **Simple and rapid bone removal** reduces the stress on the patient, thus substantially alleviating postoperative complaints.
- **Optimal positioning of transplantation material:** It is often possible to remove bone in the immediate vicinity of the implantation site.
- **The high-quality blades** enable optimum cutting of the bone material into fine discs (lamellas) without exerting a great deal of pressure.




How the best perform

Hu-Friedy

IMPLANT & MEMBRANE SURGERY KIT

Prof. Dr. Daniel Buser



FAX ORDER

Partcode	Description	Quantity
10-130-05E	Scalpel blade holder	
10-256-15	Scalpel blades #15C, 100 pcs/pkt	
PPR3/9A6	PRICHARD periosteal, handle #6	
PPBUSER6	Periosteal, handle #6	
TP5042	Adson PERMA SHARP tissue forceps surgery	
TP5041	Adson PERMA SHARP tissue forceps plain	
TP5090	Tissue pliers De Bakey, PERMA SHARP, straight, 19 cm	
CL866	Lucas bone curette, handle #6	
EXC220/16	Darby Perry excavator, english, handle #6	
SMB	Scissors Metzenbaum Baby	
H4	2 x Halsted-Mosquito artery forceps, straight	
PCV12KIT6	Colorvue Kit (7 tips, 1 handle 3-6-9-12)	
S5083	Iris Perma Sharp surgic. scissors, curved, 11,5 cm	
PH26M	Periosteal micro surgery, handle #6	
NH5052	PERMA SHARP needle holder	
MAL1	Bone mallet	
CBFL214	Bone chisel, 4 mm x 135 mm, Partsch type	
CBFL112	Bone chisel, 6 mm x 135 mm	
PPBUSE1/26	Dr. Buser scraper (3+4 mm) handle #6	
PPBUSE1/36	Dr. Buser scraper (3+5 mm) handle #6	
PPBUSE2/36	Dr. Buser scraper (4+5 mm) handle #6	
RSMID1	Middeldorpf retractor, 17 x 14 mm	
RSMID2	Middeldorpf retractor, 22 x 20 mm	
IMS-1372SH	IMS bur cushion 1/2, small lid	
IMDIN133S	IMS signature cassette for 16 instruments, orange	
PKN1526	KN 152 periosteal, handle #6	
IMP6345S6	Implant instrument no. 710 sinus lift, handle #6	
IMP65786	Kramer-Nevins sinus lift instrument, handle #6	
IMP6577SC6	Kramer-Nevins sinus lift instrument, handle #6	



In order to keep the scraper cutting edges as sharp as possible at all times, we recommend sharpening them with Arkansas Stone SS4 after every use!

SS4
Arkansas Stone



IMDIN13_S
DIN Cassette for 18 Instruments

IMDIN130S	RED	
IMDIN131S	GREY	
IMDIN132S	WHITE	
IMDIN133S	ORANGE	
IMDIN135S	YELLOW	
IMDIN136S	PINK	
IMDIN138S	BLUE	
IMDIN139S	GREEN	

Please complete and fax to Hu-Friedy: Free-Fax 00800 48 37 43 40

Name	Tel.	Fax
Street	Post Code/Town	
Course	Dental Dealer	
Signature/Stamp	Date	

Please visit us at www.Hu-Friedy.eu

How the best perform

