

CONFIDENTIAL

HU-FRIEDY® COLORVUE® PROBE
STOP SQUINTING TO READ POCKET DEPTHS!



FIRST BECAUSE WE LAST.™

Hu-Friedy Mfg. Co., Inc.
3232 N. Rockwell Street
Chicago, IL 60618
1-800-Hu-Friedy
www.hu-friedy.com

PART CODES & DESCRIPTIONS

Intro Kits

PCVPROM011	Colorvue Probe #11, Markings 3-6-8-11; 1 tip and 1 handle
PCVPROM0	Colorvue Probe #12, Markings 3-6-9-12; 1 tip and 1 handle
PCVPROMOWM	Colorvue Probe, Williams Markings 1-2-3-5-7-8-9-10; 1 tip and 1 handle
PCVPROM0NC	Colorvue Probe, UNC12 Markings 1-2-3-4-5-6-7-8-9-10-11-12; 1 tip and 1 handle

12 Packs

PCV12KIT12	Colorvue Probe #12, Markings 3-6-9-12; 12 tips and 2 handles
PCV11KIT12	Colorvue Probe #11, Markings 3-6-8-11; 12 tips and 2 handles
PCVNCKIT12	Colorvue Probe, UNC12 Markings 1-2-3-4-5-6-7-8-9-10-11-12; 12 tips and 2 handles

6 Packs

PCV12KIT6	Colorvue Probe #12, Markings 3-6-9-12; 7 tips and 1 handle
PCV11KIT6	Colorvue Probe #11, Markings 3-6-8-11; 7 tips and 1 handle
PCVWKIT6	Colorvue Probe, Williams Markings 1-2-3-5-7-8-9-10; 7 tips and 1 handle
PCVNCKIT6	Colorvue Probe, UNC12 Markings 1-2-3-4-5-6-7-8-9-10-11-12; 7 tips and 1 handle

Replacement Tips

PCV12PT	Colorvue Probe #12, Markings 3-6-9-12; 12 tips
PCV11PT	Colorvue Probe #11, Markings 3-6-8-11; 12 tips
PCVWPT	Colorvue Probe, Williams Markings 1-2-3-5-7-8-9-10; 12 tips
PCVUNC12PT	Colorvue Probe, UNC12 Markings 1-2-3-4-5-6-7-8-9-10-11-12; 12 tips

Handles


PH6	Colorvue Probe Satin Steel® Handle
XP23/PH6	#23 Colorvue Expro Handle*
XPTU17/PH6	#TU17 Colorvue Expro Handle*
XP17/PH6	#17 Colorvue Expro Handle*

*Colorvue tips sold separately

EXPRO DESIGNS

Available for clinicians who prefer the convenience of a probe and explorer combination. The Colorvue Probe is available in three popular expro designs.

*Colorvue tips sold separately



XP23/PH6
#23 Colorvue Expro*



XP17/PH6
#17 Colorvue Expro*



XPTU17/PH6
#TU17 Colorvue Expro*



XPTU17/PH6

CONTRAST

COLORVUE PROBES allow you more time for yourself and patient education!

- Vivid yellow tip & black markings provide increased intraoral visibility for faster and more accurate assessments
- Flexible, rounded tip ensures greater patient comfort and acceptability
- Convenient twist-on design allows for easy replacement of worn tips and use of larger diameter, ergonomic handle
- Safe for use around implants

Colorvue®



SEE THE VISIBLE DIFFERENCE!



1



EASIER ASSESSMENTS!



2

1. Compare the Colorvue probe to a standard metal probe. Colorvue's yellow tip with black markings provides superior contrast to the gingival tissue when measuring pocket depth.
2. 5 mm of recession is easily assessed as a result of the excellent contrast between the Colorvue probe and exposed root surface. In addition, the Colorvue probe enhances periodontal measurement of attachment loss.

Photos courtesy of Dr. James Pavlatos

TIP DESIGNS

Colorvue tips are available in four popular tip designs. Markings will wear with normal use; simply dispose of the tip. Replacement tips are available in packs of 12. Tips last approximately 30 uses.



11
3-6-8-11



12
3-6-9-12



WILLIAMS
1-2-3-5-7-8-9-10



UNC12
1-2-3-4-5-6-7-8-9-10-11-12

Colorvue Probes feature the lightweight ergonomically-designed Satin Steel® handle for the balance and control you prefer.



PH6 COLORVUE PROBE HANDLE