

# Implant/Bone and Soft Tissue Surgery Kit according to

Dr. Sascha A. Jovanovic, D.D.S., M.S.

Manufacturer: Hu-Friedy Mfg. Co., Inc.  
3232 N. Rockwell Street  
Chicago, IL 60618 · USA

EC Representative: Hu-Friedy Mfg. Co., Inc.,  
Zweigniederlassung Deutschland  
Rudolf-Diesel-Straße 8 · D-69181 Leimen  
Tel. +49 (0) 62 24 / 97 00-0  
Fax +49 (0) 62 24 / 97 00-96  
E-Mail: info@hu-friedy.de · www.hu-friedy.de



Dr. Sascha A. Jovanovic  
D.D.S., M.S.

### Periodontal Chisel Jovanovic 13K/TG




**PJO/16**  
4 mm, handle 6

**PJO/26**  
5 mm, handle 6

**PJO/36**  
6 mm, handle 6

These double-ended instruments feature straight and curved mono-beveled chisels designed to separate the periosteum, clean the bone and root surfaces and sculpt the bone.

### Sinuslift Instruments



**SKRA36**  
Sinuslift Instrument/  
Curette,  
handle 6

**IMPJO66**  
Sinuslift Instrument  
Jovanovic,  
handle 6

**PESIM16**  
Sinuslift Instrument  
Simion  
handle 6

**PFIWDS1MK**  
Tunneling Knife  
Woodson, handle 6




**PLGJO/46**  
Bone Plugger  
Jovanovic,  
handle 6

This double-ended instrument will condense a bone graft into a defect or sinus cavity.



**PPR356**  
Periosteal  
Prichard  
handle 6

May be used as a mirror as well as for tissue retraction.



**RBOE**  
Rongeurs  
Boehler,  
Double Action,  
16 cm  
Shown 65% size.

**BSJO5**  
Bone Scoop  
Jovanovic

This instrument is designed to harvest bone. Shown 75% size.




■ **Curriculum Vitae Dr. Sascha A. Jovanovic, D.D.S., M.S.**

Sascha A. Jovanovic was formally trained in Periodontics at UCLA School of Dentistry, in Implant Dentistry at Loma Linda University and in Prosthodontics at University of Aachen, Germany and holds a Master of Sciences degree in Oral Biology from UCLA.

He specializes in dental implant therapy and bone & soft tissue reconstruction, is a lecturer at UCLA Continuing Dental Education, is Professor (Hon.) at Univ. of British Columbia, division of Periodontics and Implants, Vancouver, Associate Professor (Hon.) at the University of Hong Kong Department of

Oral Surgery and is director and founder of the web based Global Institute for Dental Education at [www.globalinstituteonline.com](http://www.globalinstituteonline.com).

He is the past-President of the European Association for Osseointegration (EAO) and a past-Board member of the Osseointegration Foundation (OF). He was the program chair of the EAO 2000 annual scientific congress, recipient of the young investigator fellowship of the American Academy of Periodontology, of the Implantology award of the California Society of Periodontists and of the German Implant Society Research award.



**H4**  
Hemostat  
Halsted-  
Mosquito,  
straight, 12 cm  
Shown 85% size.

**CL866**  
Lucas Curette  
2,8 mm,  
handle 6

**CM2**  
Surgical  
Curette  
Molt  
Shown  
85%  
size.

**FX36E**  
Extraction Forcep,  
european style,  
serrated.  
Special for lower  
premolars.  
Shown 70% size.

**NH5024R**  
Needle Holder  
Perma Sharp Castroviejo,  
straight, 18 cm  
Shown 85% size.

**Scissors**

**S5080**  
Perma Sharp  
Scissors Goldman-  
Fox, straight,  
12,5 cm

Perfect for cutting  
sutures and  
membranes  
(1 blade serrated).  
Tapered beaks for  
interproximal use)  
Shown 85% size.

**S14SC**  
Super-Cut Scissors  
LaGrange,  
11,5 cm,  
double-curved  
Shown 85% size.

**CRMS**  
Retractor Universal  
Simion  
Shown 85% size.

He is an honorary member of the South African Society of Periodontology. Dr. Jovanovic's clinical research emphasizes bone reconstruction techniques and aesthetic management in dental implant therapy.


His applied research focuses on bone regeneration with different bone graft materials, rhBMP-2 and barrier membranes as well as the biology of soft tissues around teeth and implants.

He has lectured extensively worldwide and has published over 60 articles and book chapters, and one textbook titled *Color Atlas of Implantology* (Thieme Publ.).

He is on several editorial boards for scientific journals and has published a DVD-education series on advanced implant therapy. He is the current Editor in Chief of the *Implant Tribune*.

Sascha A. Jovanovic, D.D.S., M.S.  
 UCLA Continuing Dental Education  
 Global Institute for Dental Education  
 12217 Pico Blvd. · Los Angeles, CA 90064  
 Tel: (310) 820-9641 · Fax: (310) 207-9025  
 email: sascha@jovanoviconline.com

**Gracey Currettes Rigid, EverEdge Technology, handle 9**



**SG5/6R9**

Medium contra-angle for anteriors and premolars.

**SG7/8R9**

Medium contra-angle for premolars and molars, facial and lingual surfaces.

**SG11/12R9**

Angulated to reach mesial surfaces of posterior teeth.

**SG13/14R9**

Angulated to reach distal surfaces of posterior teeth.

**Mirrors and Mirror Handle**



**MHE6**  
 Mirror handle 6  
 Shown 85% size.



**M4C**  
 Mouth Mirror, #4,  
 Ø 22 mm  
 Shown 85% size.



**M5C**  
 Mouth Mirror, #5,  
 Ø 24 mm  
 Shown 85% size.



**PCPUNC156**  
 Color-Coded Probe  
 UNC 15, handle 6  
 Shown 85% size.

**Tissue Pliers**



**TP5070**  
 Tissue Pliers  
 General Perma  
 Sharp, 15 cm  
 Shown 70% size.

**TP5090**  
 Tissue Pliers  
 DeBakey  
 Perma Sharp, 15 cm  
 Shown 70% size.



**TREo60M\***  
 Trephine Bur,  
 Modified, 6 mm

Round, hollow, serrated burs for extracting bone, marked in 2 mm increments, titanium coated, adapts to contra-angle handpiece.

**Scalpel Handles and Blades**



**10-130-05E**  
 Scalpel Handle

Finely balanced, penlike handle that easily rotates and maneuvers in difficult to reach areas.



**10-256-15\***  
 Scalpel Blade 15C,  
 sterile, Box of 100



**10-255-15\***  
 Scalpel Blade 15,  
 sterile, Box of 100



**MB67\***  
 Micro Scalpel  
 Blade #67  
 12 pcs./pck



**MB69\***  
 Micro Scalpel  
 Blade #69  
 12 pcs./pck



**SHDPV**  
 Micro Scalpel  
 Handle

# Implant/Bone and Soft Tissue Surgery Kit according to

Dr. Sascha A. Jovanovic, D.D.S., M.S.

Fax Order



Manufacturer: Hu-Friedy Mfg. Co., Inc.  
3232 N. Rockwell Street  
Chicago, IL 60618 · USA

EC Representative: Hu-Friedy Mfg. Co., Inc.,  
Zweigniederlassung Deutschland  
Rudolf-Diesel-Straße 8 · D-69181 Leimen  
Tel. +49 (0) 62 24 / 97 00-0  
Fax +49 (0) 62 24 / 97 00-96  
E-Mail: info@hu-friedy.de · www.hu-friedy.de

Partcode	Description	Quantity
10-130-05E	Scalpel Handle	
10-255-15	Scalpel Blade 15 Sterile, Box of 100	
10-256-15	Scalpel Blade 15C Sterile, Box of 100	
BSJO5	Bone Scoop Jovanovic	
CL866	Surgical Curette 86 Lucas, handle 6	
CM2	Surgical Curette 2 Molt	
CRMS (Qty. 2)	Retractor Universal Simion	
FX36E	Extraction Forcep, europ. style, serrated	
H4 (Qty. 2)	Hemostat 4, Halsted-Mosquito, Straight 12cm	
IMDIN13_S	DIN cassette for 18 instruments	
IMPJO6	Implant Instrument Jovanovic	
MHE6 (Qty. 2)	Mirror handle 6	
M4C	Mouth Mirror #4, Ø 22 mm	
M5C	Mouth Mirror #5, Ø 24 mm	
MB67	Micro Scalpel Blade #67 12 pcs./pck	
MB69	Micro Scalpel Blade #69 12 pcs./pck	
NH5024R	Needle Holder Perma Sharp, Straight Castroviejo, 18cm	
PCPUNC156	Color-Coded Probe 15 UNC, handle 6	
PESIM16	Sinus Lift Instrument, handle 6	
PFIWDS1MK	Tunneling Knife Woodson, handle 6	
PJO/16	Periodontal Chisel Jovanovic 13K/TG, 4mm, handle 6	
PJO/26	Periodontal Chisel Jovanovic 13K/TG, 5mm, handle 6	
PJO/36	Periodontal Chisel Jovanovic 13K/TG, 6mm, handle 6	
PLGJO/46	Bone Plugger 4 Jovanovic, handle 6	
PPR356	Periosteal 3 Prichard, handle 6	
RBOE	Rongeurs Boehler	
S14SC	Super-Cut Scissors 14 LaGrange, 11.5cm	
S5080	Perma Sharp Scissors Goldman-Fox, Straight, 12.5cm	
SG11/12R9	Rigid Gracey Curette 11/12, handle 9, EverEdge Technology	
SG13/14R9	Rigid Gracey Curette 13/14, handle 9, EverEdge Technology	
SG5/6R9	Rigid Gracey Curette 5/6, handle 9, EverEdge Technology	
SG7/8R9	Rigid Gracey Curette 7/8, handle 9, EverEdge Technology	
SHDPV	Micro Scalpel Handle	
SKRA36	Periodontal Surgical Curette 3 Kramer-Nevins, handle 6	
TP5070	Perma Sharp Tissue Pliers General, 15cm	
TP5090	Perma Sharp Tissue Pliers DeBakey, 15cm	
TREo60M	Trephine Bur 6mm Modified	



**DIN cassette IMDIN13\_S**  
DIN cassette  
for 18 instruments  
plus additional space  
Dimensions: 290 x 184  
x 34 mm,  
with space, for surgery.

- Grey ..... IMDIN131S
- White .....IMDIN132S
- Orange.....IMDIN133S
- Yellow ..... IMDIN135S
- Pink ..... IMDIN136S
- Blue ..... IMDIN138S
- Green ..... IMDIN139S
- Red ..... IMDIN130S

Please complete



Name		
Street		
Post Code/Town		
Tel.	Fax	E-mail
Date	Signature	

Your dental supplier: