

Introducing the Micro Mini Five Gracey Curette

With ultra slender blades and increased shank rigidity, coupled with revolutionary EverEdge® technology, the Micro Mini Five Gracey Curettes will enable clinicians to precisely debride root and tooth surfaces free of all deposits, even in the deepest and most challenging periodontal pockets.

Blade is 20 % thinner than a Mini Five Gracey to further reduce tissue distention and ease subgingival insertion

Elongated terminal shank for access into deep periodontal pockets and root surfaces of 5mm and more

Increased shank rigidity compared to traditional Mini-Five Gracey curettes

Area-specific design patterned after the standard Gracey curettes offers exact adaptation and control



PART CODES:

