

sidekick™



sidekick™
Sharpening Technology

Hu-Friedy

sidekick™

The Edge of
Sharpening Technology



Hu-Friedy

Manufacturer: Hu-Friedy Mfg. Co., Inc.
3232 N. Rockwell Street · Chicago, IL 60618 · USA

EC Representative: Hu-Friedy Mfg. Co., Inc.,
Zweigniederlassung Deutschland
Rudolf-Diesel-Straße 8 · D-69181 Leimen
Tel. +49 (0) 62 24 / 97 00-0 · Fax 97 00-97
E-Mail: info@hu-friedy.de · www.hu-friedy.de

HF-CB-SID/R2-03/05/15.000

GB

Introducing sidekick™

So Simple. So Smart.


The Hu-Friedy Sidekick Sharpener makes routine maintenance sharpening of scalers and curettes fast and easy. This straightforward device guides you and your scaler to sharp, consistent results time after time.



sidekick™ The Smart Sharpener™

Reliable Results.

The Hu-Friedy Sidekick sharpener is the only power device with instrument Guide Channels and a Vertical Backstop to help control blade angulation. These "template-like" features allow virtually effortless positioning of your scalers/curettes and provide consistent sharpening results. Correct routine sharpening will help extend the life of your instruments.



1. For Scaler and Curettes Universal
2. For tips of Curettes (curve)
3. For Gracey Curettes

The Sidekick is Always Within Reach.

It's cordless! The small, compact, cordless Sidekick can be placed anywhere and requires minimal counter space – perfect for anywhere in your dental office.

See how easy and fun sharpening can be! Sharpen at your convenience with the Hu-Friedy Sidekick, The Smart Sharpener.

Sidekick Ordering Information

Sidekick Sharpening Kit: SDKKIT

Kit Includes: Sidekick Unit, User Manual & CD-ROM, Ceramic Sharpening Stone, Plastic Test Stick, Magnifying Lens and 2 AA Batteries.

Replacement Parts & Accessories

Guideplate Screw	SDKGPSCREW
Guideplate	SDKGPLATE
Ceramic Stone	SDKSTONEC
Arkansas Stone	SDKSTONEA
Sharpen-EZ™ Stone Oil	SSO
Plastic Test Sticks	PTS
Magnifying Lens	LENS



The Sidekick sharpener is designed to perform routine maintenance sharpening of scalers and curettes. It is not designed to recondition severely dull instruments. For instruments that are very dull, a professional sharpening service is recommended.

For questions about the Sidekick or Sidekick replacement parts contact your authorized Hu-Friedy dealer